

DODAVATEL:

3B instruments

lékařské nástroje

Petrovická 857
 592 31 Nové Město na Mor.
 Tel./Fax.: 0616/618 135, 0616/618 317
 E-mail: 3B@3b-instruments.cz
 Web: www.3b-instruments.cz

**46600-00**

Dermatome Handpiece

- Autoclavable
- 16.000 cpm
- 300 x 90 x 42 mm / 1170 g

**46604-00**

Spare blades
 One side cutting only.
 stainless steel

- 10 ea.
- 81 x 18,5 x 0,40 mm

**46603-00**

Charger unit for 2 battery
 (without battery)

- Mains voltage 100-240 Volts, 50-60 Hz, 0.9 A
- Charging current 1.5 A
- Charging voltage 9.6 Volts
- Protections: Short circuit / Overload / Over voltage
- Dimensions: 200 x 200 x 119 mm
- Weight: 2850 g

**46601-00**

Rechargeable Ni-MH battery

- Not autoclavable
- Battery capacity: 1.6 Ah
- Voltage: 9.6 Volts
- 136 x 32 x 32 mm / 270 g

**46602-00**


Autoclavable sterile funnel
 for aseptic battery loading

- 55 x 55 x 55 mm / 60 g

DODAVATEL:

3B instruments
lékařské nástroje

Petrovická 857
592 31 Nové Město na Mor.
Tel./Fax.:0616/618 135, 0616/618 317
E-mail: 3B@3b-instruments.cz
Web: www.3b-instruments.cz



46610-00
Skin Graft Mesher
• 205 x 165 x 120 mm / 4000 g



46611-00 Skin graft carrier > 1,5:1 (Box of 20)
46612-00 Skin graft carrier > 3:1 (Box of 20)